




# Anti-gp120/160 (SIV/mac251)

Product Name	Anti-gp120/160 (SIV/mac251)
Catalog No	IT-001-023
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	SIV gp120 protein (SIV/mac251) (Genebank accession # <a href="#">AAB00471</a> )
Applications	Western Blotting (1:200-1:1000), ELISA, Neutralization assay, etc.
Specificity	Reacts with SIV (SIV/mac251) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140 (SIV/mac251)

Note: