




Pp65 (UL83)(Cytomegalovirus)

Product Name	Pp65 (UL83)(Cytomegalovirus)
Catalog No	IT-005-020Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in 500mM NaCl, 50mM Phosphate buffer (pH 7.4), 200mM imidazole. Reconstitute the protein using 8M urea.
Protein	C-terminal 6xHis tagged Pp65 (UL83) immunodominant region (a.a. 297-510)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Store at -80°C; AVOID REPEATED FREEZE/THAW ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Pp65 (UL83) protein (a.a. 297-510)

Note: