

HA₁(H14N5)(A/mallard/Astrakhan/263/1982)



Product Name	HA ₁ (H14N5)(A/mallard/Astrakhan/263/1982)
Catalog No	IT-003-015p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (H14N5)(A/mallard/Astrakhan/263/1982) protein (a.a.18-347) (Genebank #: P26136)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H14N5)(A/mallard/Astrakhan/263/1982)