

# gp140(HIV-1/Clade AE)(Consensus)



---

Product Name	gp140(HIV-1/Clade AE) (Consensus)
Catalog No	IT-001-167p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140 (HIV-1/Clade AE)(Consensus) protein ( <a href="http://www.hiv.lanl.gov/cgi-bin/NEWALIGN/align.cgi">http://www.hiv.lanl.gov/cgi-bin/NEWALIGN/align.cgi</a> ) (a.a. 30-667). It has 98% homology to a circulating strain (Genebank #: AAG28611) (CM235) isolated in Thailand in 2000.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140 (HIV-1/Clade AE)(Consensus)