

Human TIM-3 (Ectodomain)



Product Name	Human TIM-3 (Ectodomain)
Catalog No	IT-004-TIM-3p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of human TIM-3 (Ectodomain) protein(a.a. 22-200) (Genebank #: AAL65158). TIM-3 is the member of T-cell immunoglobulin and mucin domain-containing proteins, which are the receptors for hepatitis A virus. hTIM-3 has recently been implicated as the potential cellular receptor for the Ebola virus.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human TIM-3 (Ectodomain)