

gp120(HIV-1/Clade AE)(Consensus)



Product Name	gp120(HIV-1/Clade AE)(Consensus)
Catalog No	IT-001-168p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120 (HIV-1/Clade AE)(Consensus) protein (://www.hiv.lanl.gov/cgi-bin/NEWALIGN/align.cgi)(a.a. 35-501). It has 98% homology to a circulating strain (Genebank #: AAG28611) (CM235) isolated in Thailand in 2000.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120 (HIV-1/Clade AE)(Consensus)