




Anti-gp120/160(HIV-1/Clade A)

Product Name	Anti-gp120/160 (HIV-1/Clade A)
Catalog No	IT-001-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	gp140 Env protein (a.a. 1-486) (HIV-1/Clade A/) (Genebank No.AF084657)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1(Clade A) gp120/160 proteins. Cross-reactivity to envelope protein from other clades not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing gp120 (HIV-1/Clade A)

Note: