

Anti-gp120/gp160(JRFL)(B)



Product Name	Anti- gp120/gp160(JRFL)(B)
Catalog No	IT-001-002M17
Source	Mouse monoclonal antibody, clone 21E11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120(JRFL)(B) protein (IT-001-0024p)(Genebank No. U63632)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with gp120/gp160(JRFL)(B), gp120(Con_B/2006), gp120(MG505.WOM.ENV.A2)(A), gp120(Con_A1/2006), gp120(C.ZA.1197MB)(C), gp120(DU151)(C) proteins as tested. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-002M17

