




Anti-gp120/160(HIV-1/Clade C)

Product Name	Anti- gp120/160 (Clade C)
Catalog No	IT-001-003
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 gp120(Clade C) (C.ZA.1197MB) (Genebank accession # AY463234)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with gp120/160 (Clade C). Cross-reactivity to gp120/160 from other Clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing gp120 (Clade C)

Note: