

Anti-gp120/gp140(C.ZA.1197MB)(C)



Product Name	Anti- gp120/gp140(C.ZA.1197MB)(C)
Catalog No	IT-001-003M11
Source	Mouse monoclonal antibody, clone 2D7
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing gp145 (C.ZA.1197MB)(C) protein (a.a. 24-680) (Genebank accession AY463234).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with gp120/gp140(C.ZA.1197MB)(C) as tested so far. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-003M11

