



Anti-Nef (HIV-1/Clade A)

Product Name	Anti-Nef (HIV-1/Clade A)
Catalog No	IT-001-013
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nef (HIV-1/Clade A) protein (Genebank accession #AF084657)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Nef protein from HIV-1/Clade A viruses. Cross-reactivity to Nef from other clades not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing Nef

Note: