



Anti-Nef(Pan)(HIV-1)

Product Name	Anti-Nef (Pan)(HIV-1)
Catalog No	IT-001-014PAN
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nef protein(Con-Con) (http://www.hiv.lanl.gov/cgi-bin/CONSENSUS_DOWNLOAD/ConsensusDownloader.cgi)
Applications	Western Blotting (1:200-1:1000), ELISA(1:1000-5000)
Specificity	Reacts with HIV-1 Nef protein from clade A, B, C.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

