




# Anti-p24(HIV-1/Clade A)

---

Product Name	Anti-p24 (HIV-1/Clade A)
Catalog No	IT-001-016
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 p24 matrix protein (Clade A) (Genebank No. AF084657)
Applications	Western Blotting (1:100-1:1000)
Specificity	Reacts with HIV-1(Clade A) p55 Gag protein and p24 matrix protein. Cross-reactivity to proteins from other clades not tested.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;  
E: 293 cell expressing p24(clade A)

Note: