



Anti-p24 (HIV-1/Clade A, B, C)

Product Name	Anti-p24 (HIV-1/Clade A, B, C)
Catalog No	IT-001-017M
Source	Mouse monoclonal antibody, clone 3G8, isotype IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS+40% glycerol
Immunogen	HIV-1 p24 protein (Clade B) (Genebank No. K03455)
Applications	Western Blotting, ELISA, IP, IF (Dilution needed to be determined by end user)
Specificity	Reacts with most HIV-1(Clade A, B, C) p55 Gag protein and p24 matrix protein.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

