


Anti- p24 (HIV-1/Clade C)



Product Name	Anti- p24 (HIV-1/Clade C)
Catalog No	IT-001-018
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	p24 (HIV-1/ C.ZA.1197MB) (Genebank accession AY463234)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 p24 and Gag(Clade C) proteins. Cross-reactivity to p24 and Gag from other clades not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing p24(Clade C)

Note: