




# Anti-Nef (SIV/mac239)

---

Product Name	Anti-Nef (SIV/mac239)
Catalog No	IT-001-026
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	SIV Nef protein (Genebank accession #M33262)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with SIV Nef protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;  
E: 293 cell expressing Nef (SIV/mac239)

Note: