



Anti-Rev(HIV-1/Clade B)

| | |
|----------------|--|
| Product Name | Anti-Rev(HIV-1/Clade B) |
| Catalog No | IT-001-040M |
| Source | Mouse Monoclonal Antibody |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide. |
| Immunogen | HIV-1 Rev protein (Clade B) (Genebank No. K03455) |
| Applications | Western Blotting (2-5 µg/ml) |
| Specificity | Reacts with HIV-1(Clade B) Rev protein. Cross-reactivity to Gag and p24 from other subtypes not tested. |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

