




Anti-Tat (HIV-1/Clade B)

Product Name	Anti-Tat(HIV-1/Clade B)
Catalog No	IT-001-041M
Source	Mouse Monoclonal Antibody
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Tat protein of HIV-1 (Genebank No. K03455)
Applications	Western Blotting (2-5 µg/ml)
Specificity	Reacts with HIV-1(Clade B) Tat protein. Cross-reactivity to Tat from other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing Tat (HXBc2)

Note: