




Anti-Human PD₁

Product Name	Anti-Human Program Death 1 (PD-1) Receptor
Catalog No	IT-001-051
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Human PD-1(a.a.1-167) (Genebank No. NP_005009)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human PD-1 protein. Cross-reactivity to PD-1 from other species not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing hPD-1(a.a. 1-167)

Note: Up-band may represent fully glycosylated protein