



Anti-Gag (HTLV-1)

Product Name	Anti-Gag(HTLV-1)
Catalog No	IT-001-061
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl rabbit anti-sera with 0.1% sodium Azide.
Immunogen	Gag protein of HTLV-1 virus (Genebank No. AAB20767)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Gag protein of HTLV-1 virus.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing Gag(HTLV-1) protein