



# Anti-GFP

---

Product Name	Anti-GFP
Catalog No	IT-000-010
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Green Fluorescent Protein (GFP)
Applications	Western Blotting (1:100-1:1000); Immuno-precipitation
Specificity	Reacts with Green Fluorescent Protein. Cross-reacts to most GFP derivatives
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

