


Anti-Gag (HIV-1/Clade B)



Product Name	Anti-Gag (HIV-1/Clade B)
Catalog No	IT-001-008
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 Gag(Clade B) protein (Genbank accession #K03455)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 p24 and Gag(Clade B) proteins. Cross-reactivity to p24 and Gag from other clades not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing p24(Clade B)
F: 293 cell expressing Gag(Clade B)

Note: