




Anti-CD4 (Human)

Product Name	Anti-hCD4
Catalog No	IT-001-CD4
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl rabbit anti-sera with 0.1% sodium Azide.
Immunogen	Human CD4 ectodomain protein (Genebank No. NM_000616)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and non-human primate CD4 proteins.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing human CD4.

Note: