

Anti- gp120(HIV-1)(Pan)



| | |
|----------------|--|
| Product Name | Anti- gp120(HIV-1)(Pan) |
| Catalog No | IT-001-Pan-M1 |
| Source | Mouse monoclonal antibody |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS. |
| Immunogen | Recombinant gp120 proteins from CladeA/B/C |
| Applications | ELISA, WB, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts strongly with most gp120 (Clade C), and moderately with gp120 (Clade A) and gp120(Clade B) |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |