



Anti-H₁ (H₁N₁/Pan)

Product Name	Anti-H ₁ (H ₁ N ₁ /Pan)
Catalog No	IT-003-0011
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Multiple H ₁ proteins
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human or avian H ₁ proteins from different period of history.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing.</u> The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

