




# Anti-H3(H3N2/Wis/2005)

---

Product Name	Anti-H3(A/H3N2/Wisconsin/67/X-161/2005)
Catalog No	IT-003-004WI
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HA (a.a.17-346) (A/Wisconsin/67/X-161/2005 (H3N2)) protein (Genebank accession #ABO37609)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with contemporary H3(H3N2) proteins . Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control  
E: 293 cell expressing H3(H3N2/Wis./2005)

Note: