




# Anti-HA(H<sub>5</sub>N<sub>1</sub>)(Anhui/1/05)

Product Name	Anti-HA(H5N1)(Anhui/1/05)
Catalog No	IT-003-005A
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345)(H <sub>5</sub> N <sub>1</sub> ) (A/Anhui/1/05 (H <sub>5</sub> N <sub>1</sub> ))
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H <sub>5</sub> (H <sub>5</sub> N <sub>1</sub> ). Cross-reactivity to other HAs not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;  
E: 293 cell expressing H<sub>5</sub> (Anhui/1/05) subtype;

Note: The upper band is HA, while the lower band is HA<sub>1</sub>.