




Anti-HA(H₅N₁)(Goose/Qinghai/05)

Product Name	Anti-HA(H ₅ N ₁)(Bar headed goose/ Qinghai/1A/05)
Catalog No	IT-003-005G
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345)(H ₅ N ₁)(A/ Bar-headed Goose/ Qinghai/12/05(H ₅ N ₁)) (Genebank No. ABE68927)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H ₅ (H ₅ N ₁). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing H5(Bar headed goose/
Qinghai/1A/05) subtype;

Note: The upper band is HA, while the lower band is HA₁.