



Anti-H₅(H₅N₁)

Product Name	Anti-H ₅ (H ₅ N ₁)(A/HongKong/483/97)
Catalog No	IT-003-005H
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	HA ₁ (H ₅ N ₁)(A/HongKong/483/97)(a.a.18-346) (Genebank No. AAC32099)
Applications	Western Blotting (1:200-1:2000)
Specificity	Reacts with HA(H ₅ N ₁) protein. Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing HA(H₅N₁) (A/HongKong/483/97)

Note: