

Anti-H5(A/Vietnam/1194/2004(H5N1))



Product Name	Anti- H5(A/Vietnam/1194/2004(H5N1))
Catalog No	IT-003-005M25
Source	Mouse monoclonal antibody, clone 8D3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1 (A/Vietnam/1194/2004(H5N1)) protein (a.a. 1-345) (IT-003-005p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H5 (A/Vietnam/1194/2004(H5N1)) and H5 (A/Vietnam/1203/2004(H5N1)) only. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-005M25

