



Anti-HA(H₅N₁)(VN1203/04)

Product Name	Anti-HA(H5N1)(VN1203/04)
Catalog No	IT-003-005V
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345)(H ₅ N ₁) (A/Vietnam/1203/2004(H ₅ N ₁)) (Genebank No. AY651334)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H ₅ (H5N1). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing H₅ (VN1203/04) subtype;

Note: The upper band is HA, while the lower band is HA₁.