




# Anti-HA(H5N1)(Avian)

---

Product Name	Anti-HA(H5N1)(Avian)
Catalog No	IT-003-006
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345)(H5N1)(A/ Bar-headed Goose/ Qinghai/12/05(H5N1)) (Genebank No. ABE68927)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian H5(H5N1). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control  
1: 293 cell expressing H5(1997) subtype;  
2: 293 cell expressing H5(2004) subtype;  
3: 293 cell expressing H5(2005) subtype.

Note: The upper band is HA, while the lower band is HA1.