




# Anti-HA(H6N2)

Product Name	Anti-HA(H6N2)
Catalog No	IT-003-007
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-343)(H6N2)(A/duck/Hainan/6/2004 (H6N2) ) (Genebank accession #AAV41833)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with avian H6 (H6N2). Cross-reactivity to other HAs not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control  
E: 293 cell expressing H6

Note: