

Anti-H₆(H₆N₁)



Product Name	Anti- H ₆ (A/northern shoveler/California/HKWF115/2007(H6N1))
Catalog No	IT-003-007M1
Source	Mouse monoclonal antibody, clone 85G3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing H ₄ (A/environment/Maryland/1101/2006(H4N6)) (Genebank No. ACE81692)
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H ₆ (A/northern shoveler/California/HKWF115/2007 (H6N1)) and related H ₆ subtypes. Not cross-reacts with other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

