




Anti-HA(H7N7)(Human)

Product Name	Anti-HA(H7N7)(Human)
Catalog No	IT-003-008
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hemagglutinin (amino acid 1-345) (H7N7) (A/chicken/Netherlands /1/03(H7N7))(Genebank accession #AAR02639)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian HA (H7N7). Cross-reactivity to other HAs not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing HA (H7N7)

Note: The upper band is HA, while the lower band is HA1.