




Anti-HA₂(H₅N₁)

Product Name	Anti-HA ₂ (H ₅ N ₁)
Catalog No	IT-003-010
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HA ₂ (H ₅ N ₁) protein(a.a.347-523) (Genebank accession #AY651334)
Applications	Western Blotting (1:200-1:500)
Specificity	Reacts with HA and HA ₂ (H ₅ N ₁). Cross-reactivity against other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
1: 293 cell expressing HA₂ (H₅N₁)
2: 293 cell expressing HA(H₅N₁)

Note: