



Anti-M1(H3N2)

Product Name	Anti-M1(H3N2)
Catalog No	IT-003-011
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 100 µl of rabbit anti-sera with 0.1% sodium Azide.
Immunogen	M1 protein (H3N2)(A/Wyoming/3/03) (Genebank No. AAT08000)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian M1matrix protein (H3N2). Cross-reactivity to M1 from other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

