




Anti- M1(H5N1)(Human)

Product Name	Anti- M1(H5N1)(Human)
Catalog No	IT-003-012
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Matrix protein M1 (H5N1)(A/ Vietnam/1203/2004 (H5N1)) (Genebank No. AAW80726)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian Matrix protein M1 (H5N1). Cross-reactivity to other M1 not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing M1

Note: