




Anti-H₁₂ (H₁₂N₅)

Product Name	Anti-H ₁₂ (H ₁₂ N ₅)
Catalog No	IT-003-013
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	HA ₁ of H ₁₂ N ₅ (A/pintail/Alaska/102/2005(H12N5)) (a.a 18-342) (Genebank accession #ABL67077)
Applications	Western Blotting (1:100-1:2000)
Specificity	Reacts with HA of H ₁₂ N ₅ . Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing HA(H₁₂N₅)

Note: