



Anti-H₄(H₄N₆)

Product Name	Anti-H ₄ (H ₄ N ₆)
Catalog No	IT-003-014
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	HA ₁ (H ₄ N ₆) protein(a.a.17-343) (A/environment/Maryland/1101/2006(H4N6))(Genebank accession #ABV01167)
Applications	Western Blotting (1:200-1:2000)
Specificity	Reacts with HA (H ₄ N ₆). Cross-reactivity against other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing HA(H₄N₆)

Note: