



Anti-M2 (H₁N₁)

| | |
|----------------|--|
| Product Name | Anti-M2 (H ₁ N ₁ / A/New Caledonia/20/1999(H1N1)) |
| Catalog No | IT-003-015 |
| Source | Rabbit Polyclonal Antibody |
| Formulation | Each vial contains 100 µl rabbit anti-sera with 0.1% sodium Azide. |
| Immunogen | M2 protein (H1N1) (A/New Caledonia/20/1999(H1N1)) (Genebank No. CAD30528) |
| Applications | Western Blotting (1:200-1:1000) |
| Specificity | Reacts with human and avian M2 matrix protein (H1N1). Cross-reactivity to M2 from other subtypes not tested. |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



C: 293 cell extract control;
E: 293 cell expressing M2 protein