



## Anti-H<sub>14</sub> (A/mallard/Gurjev/263/1982(H<sub>14</sub>))

---

Product Name	Anti-H <sub>14</sub> (A/mallard/Gurjev/263/1982(H <sub>14</sub> ))
Catalog No	IT-003-017
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	H <sub>14</sub> (A/mallard/Gurjev/263/1982(H <sub>14</sub> )) (Genebank accession # M35997)
Applications	Western Blotting (1:200-1:2000), ELISA, Neutralization assay, etc.
Specificity	Reacts with H <sub>14</sub> protein of influenza viruses. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H<sub>14</sub>(A/mallard/Gurjev/263/1982(H<sub>14</sub>))

Note: