



Anti-H₁₅ (A/Australia/83/00 (H₁₅N₈))

Product Name	Anti-H ₁₅ (A/Australia/83/00 (H ₁₅ N ₈))
Catalog No	IT-003-018
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	H ₁₅ (A/Australia/83/00 (H ₁₅ N ₈)) (Genebank accession # L43916).
Applications	Western Blotting (1:200-1:1000), ELISA, neutralizing assay etc.
Specificity	Reacts with H ₁₅ (H ₁₅ N ₈ /Australia/83/00) protein of influenza virus. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: H₁₅(A/Australia/83/00 (H₁₅N₈))

Note: