



Anti-NA (H₃N₂)

Product Name	Anti-NA(H ₃ N ₂)
Catalog No	IT-003-020
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NA from influenza A virus (H ₃ N ₂) (A/ Wyoming/3/03(H3N2)) (Genebank No. AAT08000)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with influenza A NA (H ₃ N ₂) protein. Cross-reactivity to other NA proteins not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

