




Anti-NA (H₅N₁)

Product Name	Anti-NA(H ₅ N ₁)
Catalog No	IT-003-021
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	NA from influenza A virus (H ₅ N ₁) (A/Vietnam/1203/2004(H5N1)) (Genebank No. AY651334)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with influenza A NA (H ₅ N ₁) protein. Cross-reactivity to other NA proteins not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C E
C: 293 cell extract control
E: 293 cell expressing NA (H₅N₁)

Note: