




Anti-NP(H1N1)

Product Name	Anti-NP(H1N1)
Catalog No	IT-003-023
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nucleocapsid protein NP (H1N1)(A/ New Caledonia/20/1999(H1N1)). (Genebank No. ABF21286)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian Nucleocapsid protein NP (H1N1). Cross-reactivity to most NP from other subtypes.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
1: 293 cell expressing NP(H1N1);
2: 293 cell expressing NP(H5N1).

Note: