



Anti-NP(Influenza A)

Product Name	Anti-NP(Influenza A)
Catalog No	IT-003-023M1
Source	Mouse monoclonal antibody, clone#65.
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG2a in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing multiple influenza A nucleoproteins.
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most NP from influenza viruses. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 3 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-023M1

