




Anti-NS1(H5N1)(Human)

Product Name	Anti- NS1(H5N1)(Human)
Catalog No	IT-003-025
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Nonstructural protein 1(NS1) (H5N1)(A/ Vietnam/1203/2004 (H5N1)) (Genebank No. AAW80934)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with human and avian Nonstructural protein 1(NS1) (H5N1). Cross-reactivity to NS1 from other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing NS1, band at *

Note: The nature of the 75kD band is unclear.