


Anti-H₁₁ (A/duck/Yangzhou/906/2002 (H₁₁N₂))



Product Name	Anti-H ₁₁ (A/duck/Yangzhou/906/2002 (H ₁₁ N ₂))
Catalog No	IT-003-035
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	H ₁₁ (A/duck/Yangzhou/906/2002 (H ₁₁ N ₂)) protein (Genebank accession AAY85533)
Applications	Western Blotting (1:200-1:1000), ELISA(1:100-1:10000)
Specificity	Reacts with H ₁₁ (H ₁₁ N ₂) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1 (A/duck/Yangzhou/906/2002 (H₁₁N₂))