

Anti-H₁₀(A/blue-winged teal/Louisiana/Sg-00073/07(H₁₀N₇))



Product Name	Anti-H ₁₀ (A/blue-winged teal/Louisiana/Sg-00073/07(H ₁₀ N ₇))
Catalog No	IT-003-034
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HA1(A/blue-winged teal/Louisiana/Sg-00073/07(H ₁₀ N ₇)) protein (a.a. 18-341) (Genebank No. ADH96304)
Applications	ELISA(1:1000-1:100,000); Neutralizing assay, IP, etc.
Specificity	Reacts with H ₁₀ (H ₁₀ N ₇). Cross-reactivity to other subtypes needs to be determined by the end user.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-034

